

World Firstly Developed Extracorporeal Kegel exercise device HnJ-350 for Home use



FURUN MEDICAL

Composition of the product



Main body



Smart phone or
Tablet PC



Low frequency
pad for
body(optional)



Low frequency
cable(Optional)



Charging cable

● Supplying package per unit


1. Main Body (Cushion) … [1EA]
2. Charging cable ……… [1EA]
3. Low Frequency pad …[4EA,Option]
4. Low frequency cable ..[1EA,Option]

Product spec

Pelvic Floor Muscles Training System (Model : HnJ-350)

1. Main Body (Cushion)	Input Power	DC 5V
	Battery Capacity	2,500mAh
	measurement range of strength	0 ~ 20kgf
	Dimensions	250 x 330 x 45mm (W-D-H)
	Weight	3kg
2. Display part	Smartphone	Android Phone / i-phone
	Tablet PC	Option :
3 Operation environment	Temperature	10°C - 40°C
	Humidity	0% - 75% RH
4. Transport and storage condition	Temperature	0°C - 40°C
	Humidity	0% - 75% RH
	Air Pressure	700 – 1060 mbar

Comparison table with global competitors

	Furun Medical		Global competitors	
	Bio feedback device with cloths on (World firstly developed)		Bio feed back device, Insertion probe into the vagina (most widely using now)	
Product type			 <p>Kegel Exercises: Contract your pelvic floor muscles for three seconds, then relax the muscles for three seconds. Do this 10-15 times several times a day. Although shown here while lying down, these exercises can be done during a variety of daily activities, such as sitting in a meeting, while stopped in your car at a traffic light or when talking on the phone.</p>	
Exercise scheme	Pelvic floor muscle training and low frequency stimulation together with cloths on		Pelvic floor muscle training with inserting probe into the vagina (Graph type)	
Convenience	Can use with cloths on / Non insertion type	0	Take off the cloths / probe insertion type	X
Motivation	Various kinds of entertainments(Games & Fun)	0	Simple graph/ number	X
Treatment way	Muscle training (Kegel exercise) Low frequency stimulation (urinary incontinence+ reduce pain)	0	Simple muscle training	△
Target	For Man, Woman, Old, Youth every body, Anywhere	0	Woman, only separated space	X

Function of the product

- Kegel exercise mode

<p>screen</p>			
	<p>Graph mode</p>	<p>Game 1</p>	<p>Game 2</p>
<p>Des cript ion</p>	<p>Graph mode appearing various short and long twitch exercise graph.</p>	<p>Designed to score coins by exercise pelvic floor muscles to improve pelvic floor muscle strength</p>	<p>Improving controlling ability for the pelvic floor muscles with contraction and relaxation.</p>

- Low-frequency stimulation mode



Low frequency screen



Pelvic stimulation



Pain Relief

Features

Same principle
with hospital
device

The Same bio feed back device principle as hospital device which are using in big hospitals, clinics now

Convenience

Realized exactly same pelvic floor muscle exercise with cloths on which was possible by inserting a probe into the female vagina in the past

Various
functions

Exercise therapy through Kegel exercise,
Nerve stimulation therapy through low frequency stimulation.
waist, shoulder, knee ,body, etc.. low frequency pain treatment

Easy to
use

Mother, father, grand mother, grand father, whoever...
Everybody can easily use.

Maximum
Treatment
effectiveness

Maximizing treatment effect by customized
exercise therapy and stimulation therapy to user.